

HPAI diagnosis and laboratory biosafety

Ji-Ming Chen, Ph.D

Head of the Laboratory of Surveillance

China Animal Health and Epidemiology Centre

2007.09.03


The views expressed in this paper/presentation are the views of the author and do not necessarily reflect the views or policies of the Asian Development Bank (ADB), or its Board of Governors, or the governments they represent. ADB does not guarantee the accuracy of the data included in this paper and accepts no responsibility for any consequence of their use. Terminology used may not necessarily be consistent with ADB official terms.

I. General Situation

China poultry annually increases by 3-8% for 3 decades

In 2005, China poultry:

- Total production: 14 billion
- Stocking population: 5 billion
- Chicken: 4.3 billion, 1/4 global



Infectious diseases emerge within this background

- Poultry production, transportation and live market increases so dramatically
- Easy for infectious diseases to spread and circulate
- Difficult to control and eliminate infectious diseases

HPAI: 01/2004 - 08/2007

- 94 HPAI outbreaks confirmed in 22 provinces of China
- 34.6 million poultries culled
- 24 human cases from 2005, mainly in southern China.

Now, well controlled in China

YEAR	Outbreaks	Died (1000)	Slaughtered (1000,000)
2004	50	129	8
2005	32	155	23
2006	10	47	3
2007 >August	2	1.2	0.06



Principles: Early, Quick, Strict

- **Early:** Discover, diagnose, report, confirm ASAP
- **Quick:** Response ASAP
- **Strict:** Strict measures

Contain the outbreak

Minimize the spread

Vaccination & Stamping-out Policy

- Vaccination: all domestic birds are compulsory for vaccination, even for backyard flocks.
- Stamping-out: all the infected and suspect animals within the epidemic spot should be slaughtered, well buried and disinfected.

II. Lab Biosafety of HPAI

Who benefits from biosafety?

- **Farmers (direct)**
- **Government (indirect)**

More tax income and less expenditure

Lab Biosafety Principle 1

Double guarantee

- Containment (Keep pathogens in)
- Disinfection (Kill pathogens)

1. Containment

Keep pathogens in

- The tips
- The tubes
- The bottles
- The boxes
- The rooms
- The buildings



2. Disinfection

- **Physical:** Boil; Autoclave; Ultraviolet;
- **Chemical:** Acid; Alkaline; Ethanol;



Lab Biosafety Principle 2

4 components:

- **Recognition**
- **Management**
- **Behaviors**
- **Facilities**

Example 1: Africa



Example 2: Hong Kong



Minimum requirements



- **Nose mask**
- **Clothes**
- **Gloves**
- **Boots**
- **Wash hands**

Two basic concepts


- Risk Classification of pathogens
- Lab Biosafety Levels

1. Risk Classification of pathogens


Risk Group	Individual risk	Community risk
1	low	low
2	moderate	low
3	high	low
4	high	high

2. Lab Biosafety Levels (BSL)


- Basic – BSL1, BSL2
- Containment – BSL3
- Maximum containment – BSL4



BSL1 LAB



BSL2 LAB



BSL3 LAB




BSL4 LAB

For more details:

Laboratory biosafety manual

Third edition



World Health Organization
Geneva
2004

Google search

Minimum requirements in the lab

- **Healthy**
- **Nose mask and protection clothes**
- **Wash hands**
- **Autoclave clinical samples and inoculated embryos before disposal**

III. Diagnosis overview

Clinical Diagnosis

- Sudden a lot of deaths
- Hemorrhage: multiple organs
- Especially: hock, brain and heart
- Difficult to differentiate from Newcastle disease virus infection



2 kinds of laboratory diagnosis


- Antibody detection: infected or vaccinated
- Pathogen detection: clinical diagnosis and outbreak confirmation

Antibody detection method

- Haemagglutination test (HA test)
- Haemagglutination inhibition test (HI test)


HA/HI test

- AIV has haemagglutination (HA) property
- Some other viruses, like NDV, also have HA property




- The HA property of **H5** subtype AIV can be inhibited by the antibodies specific to **H5** subtype AIV, but can not by the antibodies specific to **H9** subtype AIV
- HA/HI can be used to detect unknown antibodies using standard pathogens
- HA/HI can be used to detect unknown pathogens using standard antibodies


1. Place 25 μL PBS in each well




2. Make 2-fold serial dilutions of allantoic fluid



3. Each well now contains diluted allantoic fluid (except control well)



4. Add 25 μL PBS to each well




5. Add 25 μL of 1% red blood cells to each well



Leave for 45 minutes

6. Read settling patterns



The last well to show complete haemagglutination contains 1 haemagglutinating unit in 25 μL . This well contains a 1 in 64 dilution. Therefore the original allantoic fluid contained 64 haemagglutinating units in each 25 μL .

Pathogen detection methods



- Virus isolation followed with HA/HI
- RT-PCR to detect the viral nucleic acid
- Pen-side kit to detect the viral protein

1. Virus isolation

- Using SPF or AIV-unvaccinated chicken embryos
- Using throat or cloacal swabs
- Detect the virus using HA/HI test after 4-day incubation

2. RT-PCR

- Viral RNA can be amplified many times through RT-PCR
- RT-PCR specific to all subtype
- RT-PCR specific to H5 subtype



3. Pen-side kit

- Immunoassay on a filter
- Quick, simple, but not sensitive enough



For more details:

- *Manual for the laboratory diagnosis of Avian Influenza*
- chenjiming2004@yahoo.com.cn

A spiral-bound notebook with a brown cover and a white page. The spiral binding is on the left side. The text "Thank you for your attention!" is written in blue on the white page.

Thank you for your attention!